



Docket No. 55046 (70207)

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANT(S): T. C. Walsh, et al.

EXAMINER: K. M. Kerr

SERIAL NO.: 10/017,324

GROUP: 1652

FILED: December 15, 2001

FOR: METHODS FOR PREPARATION OF MACROCYCLIC MOLECULES
AND MACROCYCLIC MOLECULES PREPARED THEREBY

Mail Stop: Amendment
Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

SIR:

DECLARATION UNDER 37 CFR 1.131

The undersigned declare as follows:

1. We are co-inventors of the above-identified application assigned to the President and Fellows of Harvard College.
2. Prior to September, 2000, we had reduced to practice reactions preparing macrocyclic molecules by contacting an excised thioesterase (TE) domain with a substrate that contained a nucleophile and an activated acyl residue.
3. Prior to September, 2000, we had reduced to practice macrocyclization substrates for use in preparing macrocyclic molecules that contained a nucleophile and an activated thioester group.
4. Prior to September, 2000, such macrocyclization substrates had been contacted in an aqueous media with a purified excised TE domain under conditions conducive to macrocycle formation. As evidence thereof, attached as Exhibit 1 are selected portions of a disclosure of the subject matter of the above-identified application. The disclosure attached as Exhibit 1 was generated, and actual experimental work disclosed therein was performed, prior to September, 2000. Portions of the disclosure attached as Exhibit 1, including specific dates, have been removed.

BEST AVAILABLE COPY

Walsh, et al.
U.S.S.N. 10/017,324
Page 2

5. We heretofore further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: June 10, 2004

John W. Trauger
John W. Trauger

Date: _____

Rahul M. Kohli

Date: _____

Henning D. Mootz

Date: _____

Mohamed A. Marahiel

Date: _____

Christopher T. Walsh

Date: _____

Dirk Schwarzer

Date: _____

Michael D. Burkart



Docket No. 55046 (70207)

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANT(S): T. C. Walsh, et al.

EXAMINER: E. M. Kerr

SERIAL NO.: 10/017,324

GROUP: 1:52

FILED: December 15, 2001

FOR: METHODS FOR PREPARATION OF MACROCYCLIC MOLECULES
AND MACROCYCLIC MOLECULES PREPARED THEREBY

Mail Stop: Amendment
Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

SIR:

DECLARATION UNDER 37 CFR 1.131

The undersigned declare as follows:

1. We are co-inventors of the above-identified application assigned to the President and Fellows of Harvard College.
2. Prior to September, 2000, we had reduced to practice reactions preparing macrocyclic molecules by contacting an excised thioesterase (TE) domain with a substrate that contained a nucleophile and an activated acyl residue.
3. Prior to September, 2000, we had reduced to practice macrocyclization substrates for use in preparing macrocyclic molecules that contained a nucleophile and an activated thioester group.
4. Prior to September, 2000, such macrocyclization substrates had been contacted in an aqueous media with a purified excised TE domain under conditions conducive to macrocycle formation. As evidence thereof, attached as Exhibit 1 are selected portions of a disclosure of the subject matter of the above-identified application. The disclosure attached as Exhibit 1 was generated, and actual experimental work disclosed therein was performed, prior to September, 2000. Portions of the disclosure attached as Exhibit 1, including specific data, have been removed.

Walsh, et al.
U.S.S.N. 10/017,324
Page 2

5. We hereby further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: _____

John W. Trauger

Date: 6/29/04

Rahul M. Kohli
Rahul M. Kohli

Date: _____

Henning D. Mootz

Date: _____

Mohamed A. Marahie

Date: _____

Christopher T. Walsh

Date: _____

Dirk Schwarzer

Date: _____

Michael D. Burkart

8052_442602.1



Docket No. 55046 (70207)

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANT(S): T. C. Walsh, et al.

EXAMINER: K. M. Kerr

SERIAL NO.: 10/017,324

GROUP: 1652

FILED: December 15, 2001

FOR: METHODS FOR PREPARATION OF MACROCYCLIC
MOLECULES AND MACROCYCLIC MOLECULES
PREPARED THEREBY

Mail Stop: Amendment
Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

SIR:

DECLARATION UNDER 37 CFR 1.131

The undersigned declare as follows:

1. We are co-inventors of the above-identified application assigned to the President and Fellows of Harvard College.
2. Prior to September, 2000, we had reduced to practice reactions preparing macrocyclic molecules by contacting a excised thioesterase (TE) domain with a substrate that contained a nucleophile and an activated acyl residue.
3. Prior to September, 2000, we had reduced to practice macrocyclization substrates for use in preparing macrocyclic molecules that contained a nucleophile and an activated thioester group.
4. Prior to September, 2000, such macrocyclization substrates had been contacted in an aqueous media with a purified excised TE domain under conditions conducive to macrocycle formation. As evidence thereof, attached as Exhibit 1 are selected portions of a disclosure of the subject matter of the above-identified application. The disclosure attached as Exhibit 1 was generated, and actual experimental work disclosed therein was performed, prior to September, 2000. Portions of the disclosure attached as Exhibit 1, including specific dates, have been removed.

Walsh, et al.

U.S.S.N. 10/017,324

Page 2

5. We hereby further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: _____

John W. Trauger

Date: _____

Rahul M. KohliDate: 11.6.2004_____
H. Mootz
Henning D. MootzDate: 10.6.2004_____
M. A. Marahiel
Mohamed A. Marahiel

Date: _____

Christopher T. Walsh

Date: _____

Dirk Schwarzer

Date: _____

Michael D. Burkart

BOS2_442602.1

Walsh, et al.
U.S.S.N. 10/017,324
Page 2

5. We hereby further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: _____

John W. Trauger

Date: _____

Rahul M. Kohli

Date: _____

Henning D. Mootz

Date: _____

Mohamed A. Marahiel

Date: 6/9/04

Christopher Walsh
Christopher T. Walsh

Date: _____

Dirk Schwarzer

Date: _____

Michael D. Burkart

BOS2_442602.1

Walsh, et al.
U.S.S.N. 10/017,324
Page 2

5. We hereby further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: _____

John W. Trauger

Date: _____

Rahul M. Kohli

Date: _____

Henning D. Mootz

Date: _____

Mohamed A. Marahiel

Date: _____

Christopher T. Walsh


Date: 06/11/2004

Dirk Schwarzer
Dirk Schwarzer

Date: _____

Michael D. Burkart

BOS2_442602.1



Docket No. 55046 (70207)

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANT(S): T. C. Walsh, et al.

EXAMINER: K. M. Kerr

SERIAL NO.: 10/017,324

GROUP: 1652

FILED: December 15, 2001

FOR: METHODS FOR PREPARATION OF MACROCYCLIC MOLECULES
AND MACROCYCLIC MOLECULES PREPARED THEREBY

Mail Stop: Amendment
Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

SIR:

DECLARATION UNDER 37 CFR 1.131

The undersigned declare as follows:

1. We are co-inventors of the above-identified application assigned to the President and Fellows of Harvard College.
2. Prior to September, 2000, we had reduced to practice reactions preparing macrocyclic molecules by contacting a excised thioesterase (TE) domain with a substrate that contained a nucleophile and an activated acyl residue.
3. Prior to September, 2000, we had reduced to practice macrocyclization substrates for use in preparing macrocyclic molecules that contained a nucleophile and an activated thioester group.
4. Prior to September, 2000, such macrocyclization substrates had been contacted in an aqueous media with a purified excised TE domain under conditions conducive to macrocycle formation. As evidence thereof, attached as Exhibit 1 are selected portions of a disclosure of the subject matter of the above-identified application. The disclosure attached as Exhibit 1 was generated, and actual experimental work disclosed therein was performed, prior to September, 2000. Portions of the disclosure attached as Exhibit 1, including specific dates, have been removed.

Walsh, et al.
U.S.S.N. 10/017,324
Page 2

5. We hereby further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: _____

John W. Trauger

Date: _____

Rahul M. Kohli

Date: _____

Henning D. Mootz

Date: _____

Mohamed A. Marahiel

Date: _____

Christopher T. Walsh

Date: _____

Dirk Schwarzer

Date: 06/09/2004



Michael D. Burkart


BOS2_442602.1

Expression ① pQE6-TE in SG600 (Amp^r)
 ② pA4-H36 in BL21 (DE3) (Amp^r)
 ③ pA6-H36 in BL21 (DE3) pLysS (Amp^r/cll^r)

— overnight prealt^r 30°C
 — inoculate 15 ml/L LB + antibiotic(s)
 — shake:

Time	(1x1L) A4	(1x1L) A6	(2x1L) TE	
2 hr	.049	.22	.36	34°C → RT
2.5 hr	.058	.36	.54	26°C
				induce w/ 400 μM IPTG
3 hr	.069	.51		26°C
	flushed	induce w/ 200 μM IPTG		
7 hr, 45 min	X	2.0	2.0	26°C

— harvest, resuspend low imidazole 1x lysis buffer (2 mM imidazole)



PRE-TE purification

- from p. 16 - expressed in "3600g" cells (AraGen)
- French press (2x) Spin 15,000 rpm / 45 min.
- Shake supernatant w/ 4 mL NTA (50% shory)
- washed 1x with 40 mL buffer at 4°C for 2.5 hours
- Pellet resin; 8,000 3000 rpm / 5 min
- Transfer to small column
- Wash at 1 mL/min for 60 min. with 20 mM imidazole, 200 mM NaCl, 20 mM Tris, pH 8.
- Then elute with 20 mM → 250 mM imidazole gradient in 60 min. The product eluted at (very approx) 70 mM imidazole
- Protein-containing fractions identified by Bradford assay (4x fraction + 200x 1:4 dye)
- split into early and late (6-7 and 8-9 respectively) fractions and dialyzed overnight at 4°C vs.

50 mM Tris-HCl, pH 8
100 mM NaCl
10 mM MgCl₂
1 mM EDTA
10% glycerol

size estimation
1.7 μm ⇒ A₂₈₀ = 0.094
A = ε c l
ε = $\frac{A}{c \cdot l} = \frac{0.094}{1.7 \times 10^{-3} \text{ m} \cdot 1 \text{ cm}}$
= 55,294 & 55,000

- Quick-freeze in LN₂ and store at -80°C.
- A₂₈₀ of 1:25 dilutions:

6-7 8-9
0.094 0.073

Bradford Assay (fraction 6-7 only)		A ₅₉₅	
200x dye	0 (blank)	.490	
200x H ₂ O + SAMPLE	1x	.603	.11
	2x	.712	.21
	4x	.857	.37

conc (vs. BSA) ⇒ 1.4 μg/mL × 1000 = 1.4 mg/mL
3.1 μg/mL × 300 = 1.5 mg/mL


1.5 mg/mL

MW × 35 kD

$\frac{1.5 \text{ mg}}{1000 \text{ mL}} \times 35,000 \text{ g/mol} = 0.0004 \text{ mmol}$

$\frac{0.0004 \text{ mmol}}{35,000 \text{ g/mol}} = 1.14 \times 10^{-8} \text{ mol}$

$\frac{1.5 \text{ mg}}{1.14 \times 10^{-8} \text{ mol}} = 131,578 \text{ g/mol}$



Synthesis of tyrosine linear peptide

H₂N - D-Phe - Pro - Phe - D-Phe - Asn - Gln - Tyr - Val - Orn - Leu - COOH

Notebook file is "TLP1.NBK".

Target Peptide: length = 10, MW = 1288.521

NH₂-END-Phe-SPC-Pro-SPC-DC-Phe-SPC-Phe-SPC-Asn-SPC-Gln-SPC-Tyr-SPC-DC-Val-SPC-DC-Orn-SPC-Leu-COOH

Support substitution = 0.750 meq/g
Support quantity = 0.400 g
Excess amino acid = 3.000 x

Synthesis on 2-Cl-Trityl resin (acid-sensitive linker).

Peptide Quantity = 0.300 mmoles
Theoretical Yield = 0.387 g

0.9 mmol x 265.3 mg / mmol = 239 mg

Starting Support: FMOC-Leu-Peptide-Acid

2- HOBT only for FFP esters

Added 138 mg HOBT => 3 = HOBT activation w/ DIPC, 0.7 mmol (MW = 153.1) to each vial first.

Cycle	AA	Proto	Time	Derivative	Grams	mL	Vial
22)	SPC	N	00:07:50	System Preparation	0.000		1
21)	Orn	B3*	00:50:10	Fmoc-L-Orn(Boc)-OH ✓	Orn 0.409	3.4	2
20)	DC	H3*	00:35:10	Double Couple	0.409	3.4	3
19)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	4
18)	Val	B2	00:50:10	Fmoc-L-Val-OPfp ✓	Val 0.455	3.5	5
17)	DC	H2	00:35:10	Double Couple	0.455	3.5	6
16)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	7
15)	Tyr	B3*	00:50:10	Fmoc-L-Tyr(tBu)-OH ✓	Tyr 0.412	3.4	8
14)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	9
13)	Gln	B2	00:50:10	Fmoc-L-Gln(Trt)-OPfp ✓	Gln 0.699	3.7	10
12)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	11
11)	Asn	B2	00:50:10	Fmoc-L-Asn(Trt)-OPfp ✓	Asn 0.686	3.7	12
10)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	13
9)	Phe	B3*	00:50:10	Fmoc-D-Phe-OH ✓	D-Phe 0.349	3.3	14
8)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	15
7)	Phe	B2	00:50:10	Fmoc-L-Phe-OPfp ✓	L-Phe 0.498	3.5	16
6)	DC	H2	00:35:10	Double Couple	0.498	3.5	17
5)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	18
4)	Pro	B2	00:50:10	Fmoc-L-Pro-OPfp ✓	Pro 0.453	3.5	19
3)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	20
2)	Phe	B3*	00:50:10	Boc-D-Phe-OH 265.3	0.239 0.349	3.3	21
1)	END	G	00:07:15	Final Cycle	D-Phe 0.000		22

Minimum loop size = 10 mL
Installed loop size = 10 mL

Estimated time required for synthesis completion: 10:46:05

Estimated Reagent consumption and requirements for synthesis completion:
consumption required

Main Wash 770 mL 870 mL

Deblock	229 mL	329 mL
Wash 2	31 mL	131 mL
Aux Wash	61 mL	161 mL
Syringe 2	34 mL	46 mL
Syringe 3	24 mL	36 mL
AAM Wash	200 mL	300 mL
Synth Waste	819 mL	
SP1 Waste	56 mL	
AAM Waste	477 mL	

Page 2

Quantity

Chemistry

After synthesis, blow out resin with N_2 , dry on lyophilizer for a while (min 5 hrs, started at 11 AM).
 I will try the acetic acid cleavage method because it looks easiest (1 step). Treat resin in solvent w/ Teflon top with AcOH/TFE/DCM (2:2:6) for 2 hours at R.T. (Nava Biochem p. 553).

100 mL cleavage mix: 20 mL AcOH
 20 mL TFE
 60 mL DCM

Filter, check filtrate by TLC

Wash 3X with cleavage mix, monitor by TLC - traces (decreasing auto.) in all washes.

(Tried second 2 hour cleavage. See no additional material. \Rightarrow cleavage done in 2 hr)




Add 15 volumes hexane and rotate dry, add 50 mL hexane and dry, transfer to 25 mL

flask w/ $CHCl_3$, dry mostly, add hexane and dry to white solid. Single spot by TLC!



(10% MeOH/ $CHCl_3$) (visualize by UV). Time: 33.5757g
 Yield: 462 mg MW = (protected) 2029.5 = 0.23 mmol (76% not 100%)

Boc-D-Phe-Pro-L-Phe-D-Phe-Asn-Gly-Tyr-Val-Orn-Leu-100H
 Trt Trt Abu Boc

Thioester formation



80% TFA



(III-11B)

TILP

tyrosine linear peptide

	Ant.	mmol	Eq	MW	d
Prep 2	157 mg	0.077	1	2029.5	
N-Ac-cysteine	8.1 mg	10%L	1.1	119.2	1.12
DEC	16 mg		1	206.3	
HOBT	12 mg		1	153.1	
K ₂ CO ₃	5 mg				
THF	0.5 mL				

→ Prep 2 ←

cleared using 1:1:3 TFE/AcOH/DCM mix (~5 mL)
for 3 hours at RT. Filter, wash 2x 2 mL clearance mix.

→ Thioester formation Dry repeatedly (4x) from CH_2Cl_2 /hexane to remove AcOH
wash up; filter, ~~dry down~~ extract w/ 10% NaHCO_3 Dry high.
to remove HOBT (extract into CH_2Cl_2 or EtOAc). Then
concentrate and check the purity.

Yield: colorless (white) solid: 313 mg

we
was for
collected
purity.

0.50.93

col for

+Na⁺

51.0

4x DCC/HOBT mix: 64 mg DCC, 2 in 1 mL THF
48 mg HOBT } add 0.25 mL/run

Batch: 313 mg dissolved in 0.5 mL THF
add 0.25 mL/run

[Take 8 μL sample = 5 mg]
Sol'n becomes cloudy w/ white DCU ppt in ~5 minutes.
Washers strong. Add N-ethylcysteine after 5 min
with DCC/HOBT. Stir 3 hours at RT.

Take to dry down

Deprotect: use 80% TFA/ CH_2Cl_2 /5% N-Acetylcysteine (scavenger)
Deprotect for 2 hours at RT. Solution turned
bright yellow upon adding deprotect mix. use small
volume (~2 mL). Then add ~50 mL cold ether
dropwise, then store at 4°C for 3 hours. Filter
white powder. Recover by centrifugation, wash
3x with cold ether. Then dissolve in
20% CH_3CN /0.1% TFA and purify by prep.

HPLC (Water, monitor at 240 nm) 3 injections
efficient removal of the free acid. Phosphatase
Gradient: [to dryness. Then take up in H₂O (1 mL per
mg of peptide and phosphatase)]

Characterization: (M_n 1168)

Yield: 20 mg [M+H]⁺ colorless powder ← ~20% overall - ok

✓ MALDI-TOF MS: 1 measured exact mass = 1389.74 (calcd = 1389.7)

✓ HPLC: single peak; trace of free acid.

✓ λ_{max} : ~220 nm ($\epsilon_{\text{rel}} = 10$), ~280 nm ($\epsilon_{\text{rel}} = 1$)

looks good!

ϵ_{220} (calculated) = 1280 $\text{M}^{-1} \text{cm}^{-1}$ (ProtParam tool), (Expasy)

2nd attempt to cyclize TLPI-SMAC (II-116B)

Previously I found that Tris-HCl at pH 8 causes precipitation of TLPI-SMAC. This did not occur with HEPES-NaOH, pH 7, even after adding 140 mM NaCl. Conc of TLPI-SMAC in this trial was 7 mM.

called
"TEB9"

Step 1: Dialyze Q60-TE vs the following buffer overnight at 4°C:

enzyme (enzyme) = 40 μ M 25 mM HEPES-NaOH, pH 7.0
50 mM NaCl
10% glycerol

no ppt
labeled
"TEB9"

→ Note: I dialyzed material called "fract. 8-9" (see p. 17)

Step 2: Dissolve peptide in buffer/water (first dissolve in water, then add buffer, then NaCl: (47 mM NaCl).

1.8 mg + 140 \times H₂O

7.5 \times 0.5 M HEPES, pH 7
1.9 \times 4 M NaCl.

Step 3: Set rxns:

- "Time 0"
- ① no peptide control (enzyme only)
 - ② no enzyme control (peptide only)
 - ③ 20 \times enzyme + 20 \times peptide

Final [TLPI-SMAC] = 3.5 mM

Final [enzyme] = 20 μ M

ppt. formed, remains upon dilution with additional 60 \times H₂O. [NaCl] = 50 mM

ppt. NOT DUE TO:

- glycerol in enzyme
- TFA anion in ~~but~~ peptide prep

"Time 1" (4) To 20 λ peptide, add 60 λ H₂O, then 10 λ enzyme. Cloudy solution, clear so than before (time sample).

$$\begin{aligned}[\text{peptide}] &= 1.6 \text{ mM} \\[\text{enzyme}] &= 4 \mu\text{M} \\[\text{NaCl}] &= 17 \text{ mM}\end{aligned}$$

"Time 2" (5) 20 λ 7 mM TLA:SMAC

$$\begin{aligned}&2 \lambda \text{ 4 mM NaCl} \\&7 \lambda \text{ 0.5 M HEPES,} \\&108 \lambda \text{ H}_2\text{O} \\&3.5 \lambda \text{ TE} \\&7 \text{ ~~200~~ } \lambda\end{aligned}$$

$$\begin{aligned}[\text{peptide}] &= 1 \text{ mM} \\[\text{enzyme}] &= 1 \mu\text{M} \\[\text{NaCl}] &= 50 \text{ mM}\end{aligned}$$

slightly cloudy solution

"Time 3" (6) Same as "Time 2" but no salt added.

$$\begin{aligned}[\text{peptide}] &= 1 \text{ mM} \\[\text{enzyme}] &= 1 \mu\text{M} \\[\text{NaCl}] &= 8 \text{ mM} \\[\text{HEPES}] &= 25 \text{ mM} \\pH &= 7\end{aligned}$$

slightly cloudy ..., but clear so than Time 2 sol'n might get clearer at 37°C, -
it seemed to in ~5 min so
took all rxn out to RT

LOW SALT IS GOOD - ROOM TEMP
MAY BE PREFERABLE.

Start at noon

MAKE PEPTIDE SOLUTION IN H₂O,
allowing you to reach lower
[NaCl].

Ideal reaction conditions (I think): no NaCl,
except a little bit from the enzyme prep.

Conditions A: 10x ~~7 mM~~ 7 mM peptide in H₂O
3.5x ~~0.5 M~~ 0.5 M HEPES, pH 7
55x ~~H₂O~~ H₂O
1.0x ~~TE 89~~ TE 89
70 ~~0.75 mg~~

Peptide Stock:

0.75 mg TPI-SMAC
in 50 μ L

[HEPES] = 25 mM
[peptide-SMAC] = 1 mM
[TE] = 1 μ M
[NaCl] = 13 mM

if soluble, try to increase enzyme concentration
OR peptide-SMAC concentration.

SOLUBILITY CHECK: TPI-SMAC is readily
soluble in n-butanol \Rightarrow this is a good
solvent for extraction: extract with
n-BuOH, then concentrate to dryness. Resuspend
residue in 20% CH₃CN / 0.1% TFA/water. Run
HPLC sample. This will remove the enzyme.

n-BuOH bp = 117°C (DMF bp = 153°C)

Rollap from glass vials.

Quench ^{all} reactions (this page and pp 120-121) by
freezing at -80°C.

Isolate peptides by extraction with butanol
(2x one volume). Then concentrate
on rotavap into a glass vial using high
vacuum. Resuspend in ~~20% CH₃CN~~ 20%
CH₃CN / 0.1% TFA in H₂O.

Run HPLC. Inject 70x samples.
Gradient 0 \rightarrow 100% CH₃CN in 0.1% TFA/water.

Rxn time: ~20 hours at Room Temp

Rxn A: condition A at lpt

Rxn B: control; no enzyme

Start ~1 PM at RT.

Quench (freeze -20°C) at PM.

1 PM → 10 AM = 21 hours vs.
5 hrs (factor of 4).

Rxn C: control: no peptide

TE domain control: run set normally
w/ no peptide, diluted to 140% with
water. Inject 20% HPLC.

Rxn A: See SM covered, hydrolysis + new peak ("cycle")
(+ from TE)

Rxn B: See same hydrolysis of SM (~30%)
(no enzyme)

Rxn C: See nothing
(no peptide)

SM: Starting material purity is very good!

TE: injected TE domain on HPLC.

new peak
collected
and

submitted
for MS

New peak "cycle" on p. 124 collected off of HPLC column
MS results:

① MALDI - TOF (SEC) calc'd: 1270.65
observed: 1270.68

See data on pp 130-131.


John Frayer

John Trauger


John Trauger → 100% CH₃CN (0.1% TFA) in H₂O (0.1% TFA)
in 30 minutes,

JWT-III-124

Overlaid Traces



Overlaid Traces




Ronb50, Channel A

no arginine control

25 mM HEPES, pH 7


1 mM NACE

20 hours, RT




p-127

Roria-65, Channel A




c:\nouveau\data\john\run-65, Channel A - Time: 24.03 Min





Synthesis of peptide -SNACs

Name	Sequence	MW (TFA salt)	M+H (calc'd)	M+H (observe)
FLP-SNAC	DPhe-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1617.7	1389.7	1389.7
FLP-M1-SNAC	Phe-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1617.7	1389.7	1390.1
FLP-M2-SNAC	DAla-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1541.6	1313.7	1313.9
FLP-M3-SNAC	DPhe-Ala-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1591.7	1363.7	1363.9
FLP-M4-SNAC	DPhe-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Ala	1575.6	1461.6	
FLP-M5-SNAC	DPhe-Pro-Phe-Asn-Gln-Tyr-Val-Orn-Leu	1470.5	1242.6	
FLP-M6-SNAC	DPhe-Pro-Phe-DPhe-Asn-Ala-Gln-Tyr-Val-Orn-Leu	1688.8	1460.7	
3LP5-SNAC	DPhe-Pro-Val-Orn-Leu	918.0	690.4	690.6
3LP10-SNAC	DPhe-Pro-Val-Orn-Leu-DPhe-Pro-Val-Orn-Leu	1602.7	1260.8	

- All peptide-SNACs purified by reverse-phase C18 chromatography (HPLC) using 0.1% TFA/ acetonitrile.
 - Analyzed by MALDI-TOF MS (DCL)






Mutant + substrate Cytolization @ 5 μ M + 50 μ M.

Mutant	Am ^t	E_{280} (m ² cm ⁻¹)	MW	μ L dissolved	A_{280}	# measurements	Conc
TLP-M1	1.2 mg	1280	1617.7	148 μ L	.2331	2	4.55 mM
TLP-M2	1.0 mg	1280	1541.6	136 μ L	.2331	2	4.55 mM
TLP-M3	1.0 mg	1280	1591.7	252 μ L (2.5 mM)	.1621	3	3.17 mM
TLP-M4	2.3 mg	1280	1575.6	272 μ L	.2429	3	4.74 mM
TLP-M5	3.8 mg	1280	1470.5	317 μ L	.1962	2	3.83 mM
TLP-M6	0.7 mg	1280	1688.8	166 μ L (2.5 mM)	.1075	2	2.10 mM
GLP-GLP	2.0 mg		1602.7	200 μ L			8.73 mM

#	Name	Abs<280nm>	#	Name	Abs<280nm>
1	TLP-M1 #1	0.23282	8	TLP-M4 #1	0.23152
2	TLP-M1 #2	0.23330	9	TLP-M4 #2	0.24377
3	TLP-M2 #1	0.22852	10	TLP-M4 #3	0.25348
4	TLP-M2 #2	0.23761	11	TLP-M5 #1	0.19374
5	TLP-M3 #1	0.15737	12	TLP-M5 #2	0.19859
6	TLP-M3 #2	0.16524	13	TLP-M6 #1	0.10589
7	TLP-M3 #3	0.16367	14	TLP-M6 #2	0.10917

Spec of glp glp decapeptide 350 μ M

Overlaid Spectra:



#	Name	$E_{220} = 8646$ Abs<220nm>	Abs<240nm>	$E_{250} = 2048$ Abs<250nm>	$E_{265} = 950$ Abs<265nm>
1	GLPGLP #1	3.12700	1.62520	0.72974	0.16068
2	GLPGLP #2	3.09230	1.59180	0.71129	0.15603
3	ignore	3.07710	1.57870	0.70660	0.15548
4	GLPGLP #3	3.08120	1.57790	0.70647	0.15559

#	Name	Abs<280nm>	#	Name	Abs<280nm>
1	GLPGLP #1	3.6455E-2	3	ignore	3.5021E-2
2	GLPGLP #2	3.4674E-2	4		3.5275E-2

For cyclization rxns make up 10x stock of TLP + Mutants + GLPGLP.

	Stock	dilute stock $\rightarrow 500 \mu\text{M}$	dilute $500 \mu\text{M} \rightarrow 50 \mu\text{M}$
TLP1	3.45 mM	17.4 $\mu\text{L} \rightarrow 120 \mu\text{L}$	10 $\mu\text{L} \rightarrow 100 \mu\text{L}$
TLP-M1	4.55 mM	13.2 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"
TLP-M2	4.55 mM	13.2 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"
TLP-M3	3.17 mM	18.9 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"
TLP-M4	4.74 mM	12.7 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"
TLP-M5	3.63 mM	15.7 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"
TLP-M6	2.10 mM	28.6 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"
GLPGLP	2.73 mM	6.9 $\mu\text{L} \rightarrow 120 \mu\text{L}$	"

Make up dilution of TyCIE as per 1/21/00 to generate 80 nM stock.

Make one set for all 50 μM RXNS. and fresh for 5 μM RXNS.

RXNS:

As per 1/21 for normal rxns + 0 time points

Store samples @ -80°C.

Repeat cyclization of TLP1-SNac w/ 8 nM TyCIE.

TLP1 stock - $A_{280} = (1/2)(.22804 + .23664) = .2324$ for 1.25 dilution. = 4.54 mM

Dilutions:

dilute from 500 μM	500 μM	463 μL \rightarrow 420 μL
	350 μM	126 μL \rightarrow 180 μL
	200 μM	40 μL \rightarrow 100 μL
	150 μM	54 μL \rightarrow 180 μL
	100 μM	60 μL \rightarrow 300 μL
dilute from 100 μM	60 μM	180 μL \rightarrow 180 μL
	35 μM	35 μL \rightarrow 100 μL
	20 μM	20 μL \rightarrow 100 μL
	10 μM	10 μL \rightarrow 100 μL

Dilute TyCIE as per 1/21/00 to give 80 nM stock. Make fresh for each set of rxns.

RXNS:

Run a set with 0 time points at 60 μM 15 μM + 35 μM . Procedure as per 1/21/00.





Final results for TycC TE catalyzed peptide-SNAC cyclization experiments

Table 1. Sequences of peptide-SNAC substrates, kinetics of TycC TE-catalyzed peptide-SNAC cyclization, and exact masses of cyclic peptide products (SNAC denotes N-acetylcysteamine thioester).

Peptide-SNAC Substrate*		Cyclization			Cyclic Peptide	
		k _{cat} (min ⁻¹)	K _M (μM)	k _{cat} /K _M (μM ⁻¹ min ⁻¹)	M+H (calculated)	M+H (observed)
TLP	dPheProPheD PheAsnGlnTyrValOrnLeu-SNAC	59 ± 13	3 ± 1	21	1270.7	1270.7
dF1F	PheProPheD PheAsnGlnTyrValOrnLeu-SNAC	< 0.05†	—	—	1270.7	not detected
dF1dA	dAlaProPheD PheAsnGlnTyrValOrnLeu-SNAC	< 0.05†	—	—	1194.6	not detected
P2A	dPheAlaPheD PheAsnGlnTyrValOrnLeu-SNAC	57	3	20	1244.6	1244.4
F3A	dPheProAlaD PheAsnGlnTyrValOrnLeu-SNAC	47	6	8	1194.6	1194.7
dF4dA	dPheProPheD AlaAsnGlnTyrValOrnLeu-SNAC	105	6	16	1194.6	1194.7
N5A	dPheProPheD PheAlaGlnTyrValOrnLeu-SNAC	30	6	5	1227.7	1227.8
Q6A	dPheProPheD PheAsnAlaTyrValOrnLeu-SNAC	33	4	8	1213.6	1213.7
Y7A	dPheProPheD PheAsnGlnAlaOrnLeu-SNAC	43	15	3	1178.4	1178.8
V8A	dPheProPheD PheAsnGlnTyrAlaOrnLeu-SNAC	49	9	5	1242.6	1242.6
O9E	dPheProPheD PheAsnGlnTyrValGluLeu-SNAC	0.5	56	0.01	1285.6	1285.6
L10A	dPheProPheD PheAsnGlnTyrValOrnAla-SNAC	15	6	3	1228.6	1228.7
TLP9	dPheProPheAsnGlnTyrValOrnLeu-SNAC	4	6	0.6	1123.6	1123.8
TLP11	dPheProPheD PheAsnAlaGlnTyrValOrnLeu-SNAC	49	20	2	1341.7	1341.4
GLP10	dPheProValOrnLeuPheProValOrnLeu-SNAC	12	5	2	1141.7	1141.8

*Residues that differ from those in the wild-type substrate TLP are in bold type. †Lower limit of detection.

Many of the kinetic assays done by Rahul Kohli



Cyclization reactions: look for cyclization

Stock solutions: (3-5 mM)

3 mM ① FLPI (MW = 1484.6, TFA salt), 2.2 mg \rightarrow 3 mM
dissolve in 494 μ l H₂O \Rightarrow 3 mM

2.4 mM ② TLP3

2.43 mM stock solution previously prepared

3 mM ③ MIF (MW = 1354.4, 2 TFA salt), 3.5 mg \rightarrow 3 mM
dissolve in 862 μ l H₂O \Rightarrow 3 mM

5 mM ④ MIB (MW = 1109.1, 2 TFA salt), 3.2 mg \rightarrow 5 mM
dissolve in 962 μ l H₂O \Rightarrow 5 mM

1 mM ⑤ Phelac (MW = 504.7, 1 TFA salt), 1.4 mg \rightarrow 1 mM
dissolve in 930 μ l \Rightarrow 1 mM

Enzymes

	TE	prep	conc.
(A)	TycCTE	TE 89	40 μ M
(B)	PCP10TE	PCP1078	23 μ M
(C)	FenTE	FenTE 46	97 μ M

substrate: 2 μ M 10 μ M 50 μ M
 enzyme: 500 μ M 200 μ M TE

97

Cyclization conditions:

30 - 200 μ M peptide-SNAC
 200 μ M TE
 2.5 mM MOPS, pH 7.0
 time = 2 min. } total volume of 400 μ L

Enzymes { PCPIOTE x 1
 Fente x 1
 TyccTE x 4

Reactions

Rxn	TE	Substrate	Amt. Peptide stock	Amt H ₂ O
1	TyccTE 200 μ M	TLP3 (320)	33 \times	287 \times
2	TyccTE 200 μ M	PheIac1 (320)	80 \times	240 \times
3	TyccTE 180 μ M	M17 400	27 \times	373 \times
4	TyccTE 180 μ M	M18 400	27 \times	373 \times
5	PCPIOTE 200 μ M	TLP3 (320)	33 \times	287 \times
6	Fente 180 μ M	FLP1 400	27 \times	373 \times

Enzyme stocks:

10.6 \times 40 μ M TyccTE 5 \times 23 μ M PCPIOTE
 189 \times 180 μ M TE buff 45 \times 400 μ M TE buff
 200 \times 2 μ M TyccTE 50 \times PCPIOTE, 2 μ M


4.3 \times 92 μ M Fente TE buff = 10 mM MOPS, pH 7,
 136 \times 180 μ M TE buff 10 mM NaCl
 200 \times 2 μ M Fente

Reaction

40 \times 250 mM MOPS, pH 7
 or 320 \times H₂O + peptide-SNAC
 40 \times TE Stock
 400 \times
 2 min at RT
 - then + 25 \times 17% TFA/water quench
 - freeze LN₂
 - + 80 \times CH₃CN (17% F/w)

↓ = new products

HPLC of Samples.



Reaction of β -series w/ Tyc TE.

Stock soln of β -series peptides.

	Amt	MW	DDW	[Final]
$\beta 6$	2.1 mg	1209.4	347 μ L	5 mM
$\beta 8-1$	1.6 mg	1403.6	456 μ L	2.5 mM
$\beta 12$	2.0 mg	1827.9	438 μ L	2.5 mM
$\beta 14$	1.6 mg	2026.0	790 μ L	1 mM

\hookrightarrow purity estimated @ 50%, thus true final [P] near 500 μ M

Rxn:

Make 500 μ M stock of each peptide - SNAC.

stock =
8.72 mM

$\beta 6$ dilute 10 μ L \rightarrow 100 μ L

$\beta 8$ dilute 20 μ L \rightarrow 100 μ L

$\beta 10$ dilute 5.7 μ L \rightarrow 100 μ L

$\beta 12$ dilute 20 μ L \rightarrow 100 μ L

$\beta 14$ use stock. (purity 50% \rightarrow ~ 500 μ M).

Rxn	(500 μ M) peptide	250 mM MOPS, pH 7	DDW	+ 2 mM Tyc TE to initiate	
$\beta 6$ TE	40 μ L	40 μ L	280 μ L	40 μ L \rightarrow 1'	22.6 μ L 20 mM + 426.4 μ L DDW } 450 μ L 2 mM Tyc TE.
$\beta 6$ cont	"	"	"	-	
$\beta 8$ TE	"	"	"	40 μ L \rightarrow 1'	add 25 μ L 1.7% TFA Flash freeze
$\beta 8$ cont	"	"	"	-	
$\beta 10$ TE	"	"	"	40 μ L \rightarrow 1'	add 25 μ L 1.7% TFA FF
$\beta 10$ cont	"	"	"	-	
$\beta 12$ TE	"	"	"	40 μ L \rightarrow 1'	add 25 μ L 1.7% TFA FF
$\beta 12$ cont	"	"	"	-	
$\beta 14$ TE	"	"	"	40 μ L \rightarrow 1'	add 25 μ L 1.7% TFA FF
$\beta 14$ cont	"	"	"	-	

For all cont rxn. add first 25 μ L 1.7% TFA.

then add 40 μ L 2 mM Tyc TE
flash freeze instantly.

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.