



# Best Available Copy

Docket No. 55046 (70207)

## IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANT(S): T. C. Walsh, et al.

EXAMINER: K. M. Kerr

SERIAL NO.: 10/017,324

GROUP: 1652

FILED: December 15, 2001

FOR: METHODS FOR PREPARATION OF MACROCYCLIC MOLECULES  
AND MACROCYCLIC MOLECULES PREPARED THEREBY

Mail Stop: Amendment  
Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

SIR:

### DECLARATION UNDER 37 CFR 1.131

The undersigned declare as follows:

1. We are co-inventors of the above-identified application assigned to the President and Fellows of Harvard College.
2. Prior to September, 2000, we had reduced to practice reactions preparing macrocyclic molecules by contacting a excised thioesterase (TE) domain with a substrate that contained a nucleophile and an activated acyl residue.
3. Prior to September, 2000, we had reduced to practice macrocyclization substrates for use in preparing macrocyclic molecules that contained a nucleophile and an activated thioester group.
4. Prior to September, 2000, such macrocyclization substrates had been contacted in an aqueous media with a purified excised TE domain under conditions conducive to macrocycle formation. As evidence thereof, attached as Exhibit 1 are selected portions of a disclosure of the subject matter of the above-identified application. The disclosure attached as Exhibit 1 was generated, and actual experimental work disclosed therein was performed, prior to September, 2000. Portions of the disclosure attached as Exhibit 1, including specific dates, have been removed.

Walsh, et al.  
U.S.S.N. 10/017,324  
Page 2

5. We hereby further declare that all statements made herein are of our own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, and that such willful false statements may jeopardize the validity of the above-identified application or any patent issued thereon.

Date: \_\_\_\_\_

John W. Trauger

Date: \_\_\_\_\_

Rahul M. Kohli

Date: \_\_\_\_\_

Hemming D. Mootz

Date: \_\_\_\_\_

Mohamed A. Marahiel

Date: \_\_\_\_\_

Christopher T. Walsh

Date: \_\_\_\_\_

Dirk Schwarzer

Date: 06/09/2004

Michael D. Burkart

BOS2\_442602.1



- Expression
- ① pQE6-TE in SG600 (Amp<sup>r</sup>)
  - ② pA4-HIS6 in BL21 (DE3) (Amp<sup>r</sup>)
  - ③ pA6-HIS6 in BL21 (DE3) pLysS (Amp<sup>r</sup>/cll<sup>r</sup>)

— overnight prealtu 30°C  
 — inoculate 15 mL/L LB + antibiotic(s)  
 — shake:

Time	(1x1L) A4	(1x1L) A6	(2x1L) TE	
2 hr	.049	.22	.36	34°C → RT
2.5 hr	.058	.36	.54	26°C


induce w/  
400λ 1m IPTG

3 hr. .069 (circled .51) 26°C

frushed induce w/  
200λ 1m IPTG

7 hr, 45 min X 2.0 2.0 26°C

— harvest, resuspend low imidazole 1x lysS buffer (2 mM imidazole)



PROTEIN PURIFICATION


- from p. 16 - expressed in "36<sup>00</sup>g (PREP4) cells (Qiagen)
- French press (2x) spin 15,000 rpm / 45 min.
- Shake supernatant w/ 4 mL N-NTA (50% storage)
- washed 1x with 40 mL buffer at 4°C for 2.5 hours
- Pellet resin; 8,000 rpm / 5 min
- Transfer to small column
- Wash at 1 mL/min for 60 min. with 20 mM imidazole, 200 mM NaCl, 20 mM Tris, pH 8.
- Then elute with 20 mM → 250 mM imidazole gradient in 60 min. The product eluted at (very approx) 70 mM imidazole
- Protein - containing fractions identified by Bradford assay (4 fractions + 200x 1:4 dye)
- split into early and late (6-7 and 8-9 respectively) fractions and dialyzed overnight at 4°C vs.

50 mM Tris-HCl, pH 8  
 100 mM NaCl  
 10 mM MgCl<sub>2</sub>  
 1 mM EDTA  
 10% glycerol

← 280 estimation  
 $1.7 \mu\text{M} \Rightarrow A_{280} = 0.094$   
 $A = \epsilon \cdot c \cdot l$   
 $\epsilon = \frac{A}{c \cdot l} = \frac{0.094}{1.7 \times 10^{-3} \text{ m} \cdot 1 \text{ cm}}$   
 $= 55,294 \times 10^3 \text{ M}^{-1} \text{ cm}^{-1}$

- Quick freeze in LN<sub>2</sub> and store at -80°C.
- A<sub>280</sub> of 1:25 dilutions:

6-7      8-9  
 0.094      0.073



Bradford Assay (fraction 6-7 only)


	Blank	A <sub>595</sub>
200x dye	0.490	
200x H <sub>2</sub> O + SAMPLE	0.603	.11
	0.712	.21
	0.857	.37

conc (vs BSA)  $\approx 1.4 \mu\text{g/mL} \times 1000 = 1.4 \text{ mg/mL}$   
 $3.1 \mu\text{g/mL} \times 1000 = 3.1 \text{ mg/mL}$

1.5 mg/mL

mw  $\approx 35 \text{ kD}$

$\frac{1.5 \text{ mg}}{1000 \text{ mL}} \times \frac{35,000 \text{ g/mol}}{1000000} = 0.04 \text{ mM}$   
 $\frac{3.1 \text{ mg}}{1000 \text{ mL}} \times \frac{35,000 \text{ g/mol}}{1000000} = 0.11 \text{ mM}$   
 $\frac{0.04 \text{ mM}}{0.11 \text{ mM}} = 0.36$



# Synthesis of tyrosine linear peptide

H<sub>2</sub>N - D-Phe - Pro - Phe - D-Phe - Asn - Gln - Tyr - Val - Orn - Leu - COOH

Notebook file is "TLP1.NBK".

Target Peptide: length = 10, MW = 1288.521

NH<sub>2</sub>-END-Phe-SPC-Pro-SPC-DC-Phe-SPC-Phe-SPC-Asn-SPC-Gln-SPC-Tyr-SPC-DC-Val-SPC-DC-Orn-SPC-Leu-COOH

Support substitution = 0.750 meq/g  
Support quantity = 0.400 g  
Excess amino acid = 3.000 x

} Synthesis on 2-Cl-Trityl resin  
(acid-sensitive linker).

Peptide Quantity = 0.300 mmoles  
Theoretical Yield = 0.387 g

0.9 mmole x 265.3 mg/mmol = 239 mg

Starting Support: FMOC-Leu-Peptide-Acid

Added 138 mg HOBT  $\Rightarrow$  2 = HOBT only for FFP esters  
3 = HOBT activation w/ DIPC<sub>1</sub>  
0.9 mmole (MW = 153.1) to each vial first.

Cycle	AA	Proto	Time	Derivative	Grams	mL	Vial
22)	SPC	N	00:07:50	System Preparation	0.000		1
21)	Orn	B3*	00:50:10	Fmoc-L-Orn(Boc)-OH ✓	Orn 0.409	3.4	2
20)	DC	H3*	00:35:10	Double Couple	0.409	3.4	3
19)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	4
18)	Val	B2	00:50:10	Fmoc-L-Val-OPfp ✓	Val 0.455	3.5	5
17)	DC	H2	00:35:10	Double Couple	0.455	3.5	6
16)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	7
15)	Tyr	B3*	00:50:10	Fmoc-L-Tyr(tBu)-OH ✓	Tyr 0.412	3.4	8
14)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	9
13)	Gln	B2	00:50:10	Fmoc-L-Gln(Trt)-OPfp ✓	Gln 0.699	3.7	10
12)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	11
11)	Asn	B2	00:50:10	Fmoc-L-Asn(Trt)-OPfp ✓	Asn 0.686	3.7	12
10)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	13
9)	Phe	B3*	00:50:10	Fmoc-D-Phe-OH ✓	D-Phe 0.349	3.3	14
8)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	15
7)	Phe	B2	00:50:10	Fmoc-L-Phe-OPfp ✓	L-Phe 0.498	3.5	16
6)	DC	H2	00:35:10	Double Couple	0.498	3.5	17
5)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	18
4)	Pro	B2	00:50:10	Fmoc-L-Pro-OPfp ✓	Pro 0.453	3.5	19
3)	SPC	J	00:09:15	N-Acetylimidazole	0.099	3.1	20
2)	Phe	B3*	00:50:10	Boc-D-Phe-OH <sup>265.3</sup>	0.239	3.3	21
1)	END	G	00:07:15	Final Cycle	0.000		22

Minimum loop size = 10 mL

Installed loop size = 10 mL

Estimated time required for synthesis completion: 10:46:05

Estimated Reagent consumption and requirements for synthesis completion:  
consumption required

MF Main Wash 770 mL 870 mL

Deblock	229 mL	329 mL
Wash 2	31 mL	131 mL
Aux Wash	61 mL	161 mL
Syringe 2	34 mL	46 mL
Syringe 3	24 mL	36 mL
AAM Wash	200 mL	300 mL
Synth Waste	819 mL	
SP1 Waste	56 mL	
AAM Waste	477 mL	

Page 2

Quantity

Chemical

- After synthesis, blow out resin with  $N_2$ , dry on lyophilizer for a while (with 5 hrs started at 11 AM).  
 - I will try the acetic acid cleavage method because it looks easiest (1 step). Treat resin in solvent w/ Teflon top with AcOH/TFE/DCM (2:2:6) for 2 hours at R.T. (Nova Biochem p. 555).




100 mL cleavage mix: 20 mL AcOH  
 20 mL TFE  
 60 mL DCM

Filter, check filtrate by TLC  
 Wash 3X with cleavage mix, monitor by TLC -  
 traces (decreasing amounts) in all washes.  
 (Tried second 2 hour cleavage, see no additional material  $\Rightarrow$  cleavage done in 2 hr)

Add 15 volumes hexane and rotate dry, add 50 mL hexane and dry, Transfer to 25 mL flask w/ chills, dry mostly, add hexane and dry to white solid. Single spot by TLC.  
 (10% MeOH/CHCl<sub>3</sub>) (visualize by UV). Tare: 33.5757g  
 Yield: 462 mg MW = (protected) 2029.5 = 0.23 mmol (76% not 60%)  
 Boc-D-Phe-Pro-L-Phe-D-Phe-Asn-Gln-Tyr-Val-Orn-Leu-COOH  
 Tst Tst Abu Boc

Thioester formation

(prep 2)



TLP  
tyrosine linear peptide

	Ant.	mmol	Eq	MW	d
Prep 2	157 mg	0.077	1	2022.5	
NAC-glycine	8.1 10%L		1.1	119.2	1.12
DCC	16 mg		1	206.3	
HOBT	12 mg		1	153.1	
K <sub>2</sub> CO <sub>3</sub>	5 mg				
THF	0.5 mL				

Prep 2

cleared mix 1:1:3 TFE/AcOH/DCM mix (~5 mL)  
for 3 hours at RT. Filter, wash 2x 2 mL cleared mix.

THF ester formation Dry repeatedly (4x) from  $\text{CH}_2\text{Cl}_2$ /hexane to remove AcOH  
wash up; filter, ~~extract~~ <sup>dry down</sup> extract w/ 10%  $\text{NaHCO}_3$  to remove HOBT (extract into  $\text{CH}_2\text{Cl}_2$  or EtOAc). Then  
concentrate and check the purity.

Yield: colorless (white) solid: 313 mg

4x DCC/HOBT mix: 64 mg DCC, 7 min 1 mL THF  
48 mg HOBT } add 0.25 mL/run

313 mg dissolved in 0.5 mL THF  
add 0.25 mL/run

Take 8  $\mu\text{L}$  sample = 5 mg  
Sol'n became cloudy w/ white DCU ppt in 2.5 minutes.  
W/ 10%  $\text{NaHCO}_3$ . Add NAC to cysteine after 5 min  
with DCC/HOBT. Stir 3 hours at RT.

Deprotect: use 80% TFA/ $\text{CH}_2\text{Cl}_2$ /5% N Acetylcysteine (scavenger)  
Deprotect for 2 hours at RT. Solution turned  
bright yellow upon adding deprotect mix. use small  
volume (~2 mL). Then add ~50 mL cold ether  
dropwise, then store at 4°C for 3 hours. Fine  
white powder. Recover by centrifugation, wash  
3x with cold ether. Then dissolve in  
20%  $\text{CH}_3\text{CN}$ /0.1% TFA and purify by prep.

HPLC (Water, monitor at 240 nm) 3 injections  
efficient removal of free free acid. Lyophilize  
Gradient: to dryness. then take up in H<sub>2</sub>O (1 mL per  
mg of peptide and lyophilize).

Characterization: (M+H)<sup>+</sup> 1168

Yield: 20 mg [M+H]<sup>+</sup> colorless powder ← ~20% overall - ok

MALDI-TOF MS: measured exact mass = 1389.74 (calcd = 1389.7)

HPLC: single peak; trace of free acid.

$\lambda_{\text{max}}$ : ~220 nm ( $\epsilon_{\text{rel}} = 10$ ), ~280 nm ( $\epsilon_{\text{rel}} = 1$ )

looks good!

$\epsilon_{220}$  (calculated) =  $1280 \text{ M}^{-1} \text{ cm}^{-1}$  (ProtParam tool), (Expasy)



## 2nd attempt to cyclize TLPI-SMAC (M-116B)

Previously I found that Tris-HCl at pH 8 causes precipitation of TLPI-SMAC. This did not occur with HEPES-NaOH, pH 7, even after adding 140 mM NaCl. Conc of TLPI-SMAC in this trial was 7 mM.

called  
"TEB"

Step 1: Dialyze Q60-TE vs the following buffer overnight at 4°C:

osmotic (enzyme) = 40 mM      25 mM HEPES-NaOH, pH 7.0  
50 mM NaCl  
10% glycerol

was ppts  
labeled  
"TEB"

→ Note: I dialyzed material called "fract. 8-9" (see p. 7)

Step 2: Dissolve peptide in buffer/water (first dissolve in water, then add buffer, then NaCl: (47 mM NaCl),

1.8 mg + 140  $\times$  H<sub>2</sub>O

7.5  $\times$  0.5 M HEPES, pH 7  
1.9  $\times$  4 M NaCl.

Step 3: Set rxns:

- ① no peptide control (enzyme only)
- ② no enzyme control (peptide only)
- "Time 0" ③ 20  $\times$  enzyme + 20  $\times$  peptide

Final [TLPI-SMAC] = 3.5 mM

Final [enzyme] = 20  $\mu$ M

ppt. forms, seems upon dilution with additional 60  $\times$  H<sub>2</sub>O. [NaCl] = 50 mM

ppt. NOT DUE TO:

- glycerol in enzyme

- TFA anion in ~~left~~ peptide prep

"Time 1" ④ To 20  $\lambda$  peptide, add 60  $\lambda$  H<sub>2</sub>O, then  
10  $\lambda$  enzyme. Cloudy solution, less so than  
before (Time 0 sample).

$$\begin{aligned}[peptide] &= 1.6 \text{ mM} \\ [enzyme] &= 4 \mu\text{M} \\ [NaCl] &= 17 \text{ mM}\end{aligned}$$

"Time 2" ⑤ 20  $\lambda$  7 mM TLP:SMAC  
2  $\lambda$  4  $\mu$  M NaCl  
7  $\lambda$  0.5 M HEPES,  
108  $\lambda$  H<sub>2</sub>O  
3.5  $\lambda$  TE  
740  $\lambda$

$$\begin{aligned}[peptide] &= 1 \text{ mM} \\ [enzyme] &= 1 \mu\text{M} \\ [NaCl] &= 50 \text{ mM}\end{aligned}$$

slightly cloudy solution

"Time 3" ⑥ Same as "Time 2" but no salt added.

$$\begin{aligned}[peptide] &= 1 \text{ mM} \\ [enzyme] &= 1 \mu\text{M} \\ [NaCl] &= 8 \text{ mM} \\ [\text{HEPES}] &= 2.5 \text{ mM} \\ \text{pH} &= 7\end{aligned}$$

slightly cloudy... but less so than Time 2  
sol'n might get clearer at 37°C, -  
it seemed to in ~5 min so  
took all rxn out to RT.

[ LOW SALT IS GOOD - ROOM TEMP  
MAY BE PREFERABLE.  
Start at noon

MAKE PEPTIDE SOLUTION IN H<sub>2</sub>O,  
allowing you to reach lower  
[NaCl].

Ideal reaction conditions (I think): no NaCl,  
except a little bit from the enzyme prep.

Conditions A: 10x ~~10x~~ 7 mM peptide in H<sub>2</sub>O  
3.5x ~~10x~~ 0.5 M HEPES, pH 7  
55x ~~10x~~ H<sub>2</sub>O  
1.0x ~~10x~~ TE89  
70 ~~10x~~

Peptide Stock:

0.75 mg  
50  $\mu$ L

TEPI-SMAC

[HEPES] = 25 mM

[peptide-SMAC] = 1 mM

[TE] = 1  $\mu$ M

[NaCl] = 1.3 mM

if soluble, try to increase enzyme concentration  
OR peptide-SMAC concentration.

SECURITY CHECK: TEPI-SMAC is readily  
soluble in n-butanol  $\Rightarrow$  this is a good  
solvent for extraction: extract with  
n-BuOH, then concentrate to dryness. Resuspend  
residue in 20% CH<sub>3</sub>CN / 0.1% TFA/water. Run  
HPLC sample. This will remove the enzyme.

n-BuOH bp = 120°C (DMF bp = 153°C)

Rotavap from glass vials.

Overexposed reactions (this page and pp 120-121) by  
freezing at -80°C.

Isolate peptide by extraction with butanol  
(2x one volume). Then concentrate  
on rotavap into a glass vial using high  
vacuum. Resuspend in ~~20%~~ 20% ~~CH<sub>3</sub>CN~~  
CH<sub>3</sub>CN / 0.1% TFA in H<sub>2</sub>O.

Run HPLC. Inject 70x samples.  
Gradient 0  $\rightarrow$  100% CH<sub>3</sub>CN in 0.1% TFA/water.

Rxn time: ~20 hours at Room Temp

Rxn A: condition A at lgt

Rxn B: control; no enzyme

Start ~1 PM at RT.  
Quench (freeze -20°C) at PM.

1 PM → 10 AM = 21 hours vs.  
5 hrs (factor of 4).

Rxn C: control no peptide

TE domain control: run set normally  
w/ no peptide, diluted to 140% with  
water. Inject 20% HPLC.

Rxn A: See SM consumed, hydrolysis + new peak ("cycle")  
(+ from TE)

Rxn B: See same hydrolysis of SM (~30%)  
(no enzyme)

Rxn C: See nothing  
(no peptide)

SM: starting material purity is very good!

TE: injected TE domain on HPLC.

new peak  
collected  
and

submitted  
for MS

New peak "cycle" on p. 124 collected off of HPLC column  
MS results:

① MALDI-TOF (SEC) calc'd: 1270.65  
observed: 1270.68

See data on pp 120-131.

John Trauger


1001

John Trauger


John Trauger → 100% CH<sub>3</sub>CN (0.1% TFA) in H<sub>2</sub>O (0.1% TFA)  
in 30 minutes,

JWT-III-124


Overlaid Traces




Overlaid Traces



R01b50, Channel A



Roria-65, Channel A






128

p. 128


c:\nouveau\data\john\loma-65, Channel A - Time: 24.03 Min





# Synthesis of peptide-SNACs

Name	Sequence	MW (TFA salt)	M+H (calc'd)	M+H (observe)
PLP-SNAC	DPhe-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1617.7	1389.7	1389.7
PLP-M1-SNAC	Phe-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1617.7	1389.7	1390.1
PLP-M2-SNAC	DAla-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1541.6	1313.7	1313.9
PLP-M3-SNAC	DPhe-Ala-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Leu	1591.7	1363.7	1363.9
PLP-M4-SNAC	DPhe-Pro-Phe-DPhe-Asn-Gln-Tyr-Val-Orn-Ala	1575.6	1461.6	
PLP-M5-SNAC	DPhe-Pro-Phe-Asn-Gln-Tyr-Val-Orn-Leu	1470.5	1242.6	
PLP-M6-SNAC	DPhe-Pro-Phe-DPhe-Asn-Ala-Gln-Tyr-Val-Orn-Leu	1688.8	1460.7	
3LP5-SNAC	DPhe-Pro-Val-Orn-Leu	918.0	690.4	690.6
3LP10-SNAC	DPhe-Pro-Val-Orn-Leu-DPhe-Pro-Val-Orn-Leu	1602.7	1260.8	

- All peptide-SNACs purified by reverse-phase C18 chromatography (HPLC) using 0.1% TFA/ acetonitrile.  
 - Analyzed by MALDI-TOF MS (DECI)






Mutant + substrate Cyclization @ 5 $\mu$ M + 50 $\mu$ M.

Mutant	Am <sup>t</sup>	$\epsilon_{280}$ (m <sup>2</sup> cm <sup>-1</sup> )	MW	$\mu$ L dissolved	A <sub>280</sub>	# measurements	Conc
TLP-M1	1.2 mg	1280	1617.7	148 $\mu$ L	.2331	2	4.55 $\mu$ M
TLP-M2	1.0 mg	1280	1541.6	130 $\mu$ L	.2331	2	4.55 $\mu$ M
TLP-M3	1.0 mg	1280	1591.7	252 $\mu$ L (25 $\mu$ M)	.1621	3	3.17 $\mu$ M
TLP-M4	2.3 mg	1280	1575.6	272 $\mu$ L	.2429	3	4.74 $\mu$ M
TLP-M5	3.8 mg	1280	1470.5	317 $\mu$ L	.1962	2	3.83 $\mu$ M
TLP-M6	0.7 mg	1280	1688.8	166 $\mu$ L (25 $\mu$ M)	.1075	2	2.10 $\mu$ M
GLP-GLP	2.8 mg		1602.7	200 $\mu$ L			8.73 $\mu$ M

#	Name	Abs<280nm>	#	Name	Abs<280nm>
1	TLPM1 #1	0.23282	8	TLPM4 #1	0.23152
2	TLPM1 #2	0.23330	9	TLPM4 #2	0.24377
3	TLPM2 #1	0.22852	10	TLPM4 #3	0.25348
4	TLPM2 #2	0.23761	11	TLPM5 #1	0.19374
5	TLPM3 #1	0.15737	12	TLPM5 #2	0.19859
6	TLPM3 #2	0.16524	13	TLPM6 #1	0.10589
7	TLPM3 #3	0.16367	14	TLPM6 #2	0.10917

Spec of glpgle decapeptide 350  $\mu$ M

## Overlaid Spectra:



1:25 dilution	For 1 M:	$\epsilon_{220} = \frac{8846}{2046}$	$\epsilon_{250} = \frac{450}{2046}$	$\epsilon_{265} = \frac{450}{2046}$	
#	Name	Abs<220nm>	Abs<240nm>	Abs<250nm>	Abs<265nm>
1	GLPGLP #1	3.12700	1.62520	0.72974	0.16068
2	GLPGLP #2	3.09230	1.59180	0.71129	0.15603
3	ignore	3.07710	1.57870	0.70660	0.15546
4	GLPGLP #3	3.08120	1.57790	0.70647	0.15559
#	Name	Abs<280nm>	#	Name	Abs<280nm>
1	GLPGLP #1	3.6455E-2	3	ignore	3.5021E-2
2	GLPGLP #2	3.4674E-2	4		3.5275E-2

For cyclization rxns make up 10x stock of TLP + Mutants + GLPGLP.

	Stock	dilute stock $\rightarrow$ 500 $\mu$ M	dilute 500 $\mu$ M $\rightarrow$ 50 $\mu$ M
TLP1	3.45 mM	17.4 $\mu$ L $\rightarrow$ 120 $\mu$ L	10 $\mu$ L $\rightarrow$ 100 $\mu$ L
TLP-M1	4.55 mM	13.2 $\mu$ L $\rightarrow$ 120 $\mu$ L	"
TLP-M2	4.55 mM	13.2 $\mu$ L $\rightarrow$ 120 $\mu$ L	"
TLP-M3	3.17 mM	18.9 $\mu$ L $\rightarrow$ 120 $\mu$ L	"
TLP-M4	4.74 mM	12.7 $\mu$ L $\rightarrow$ 120 $\mu$ L	"
TLP-M5	3.83 mM	15.8 $\mu$ L $\rightarrow$ 120 $\mu$ L	"
TLP-M6	2.10 mM	28.6 $\mu$ L $\rightarrow$ 120 $\mu$ L	"
GLPGLP	8.73 mM	6.9 $\mu$ L $\rightarrow$ 120 $\mu$ L	"

Make up dilution of TyCE as per 1/2100 to generate 80 nM stock.

Make one set for all 50  $\mu$ M RXNS and fresh for 5  $\mu$ M RXNS.

RXNS:

As per 1/21 for normal rxns + 0 time points

Store samples @ -80°C.

Repeat cyclization of TLP1-SNac w/ 8 nM TyCE.

TLP1 stock -  $A_{280} = (1/2)(.2280 + .2366) = .2324$  for 1% solution = 4.54 mM

Dilutions:

	500 $\mu$ M	463 $\mu$ L $\rightarrow$ 420 $\mu$ L
	350 $\mu$ M	126 $\mu$ L $\rightarrow$ 180 $\mu$ L
dilute from 500 $\mu$ M	200 $\mu$ M	40 $\mu$ L $\rightarrow$ 100 $\mu$ L
	150 $\mu$ M	54 $\mu$ L $\rightarrow$ 180 $\mu$ L
	100 $\mu$ M	60 $\mu$ L $\rightarrow$ 200 $\mu$ L
	20 $\mu$ M	100 $\mu$ L $\rightarrow$ 180 $\mu$ L
dilute from 100 $\mu$ M	35 $\mu$ M	35 $\mu$ L $\rightarrow$ 100 $\mu$ L
	20 $\mu$ M	20 $\mu$ L $\rightarrow$ 100 $\mu$ L
	10 $\mu$ M	10 $\mu$ L $\rightarrow$ 100 $\mu$ L

Dilute TyCE as per 1/2100 to give 80 nM stock. Make fresh for each set of rxns.

RXNS:

Run a set with 0 time points at 5  $\mu$ M  
15  $\mu$ M + 35  $\mu$ M. Procedure as per 12.102.

# Final results for TyC TE catalyzed peptide-SNAC cyclization experiments





Table 1. Sequences of peptide-SNAC substrates, kinetics of TyC TE-catalyzed peptide-SNAC cyclization, and exact masses of cyclic peptide products (SNAC denotes N-acetylcysteamine thioester).

MALDI-TOF MS

Peptide-SNAC Substrate*	Cyclization			Cyclic Peptide	
	$k_{cat}$ ( $\text{min}^{-1}$ )	$K_M$ ( $\mu\text{M}$ )	$k_{cat}/K_M$ ( $\mu\text{M}^{-1}\text{min}^{-1}$ )	M+H (calculated)	M+H (observed)
TLP	$59 \pm 13$	$3 \pm 1$	21	1270.7	1270.7
<b>dF1F</b>	$< 0.05^\dagger$	—	—	1270.7	not detected
<b>dF1dA</b>	$< 0.05^\dagger$	—	—	1194.6	not detected
P2A	57	3	20	1244.6	1244.4
F3A	47	6	8	1194.6	1194.7
<b>dF4dA</b>	105	6	16	1194.6	1194.7
N5A	30	6	5	1227.7	1227.8
Q6A	33	4	8	1213.6	1213.7
Y7A	43	15	3	1178.4	1178.8
V8A	49	9	5	1242.6	1242.6
O9E	0.5	56	0.01	1285.6	1285.6
L10A	15	6	3	1228.6	1228.7
TLP9	4	6	0.6	1123.6	1123.8
TLP11	49	20	2	1341.7	1341.4
GLP10	12	5	2	1141.7	1141.8

\*Residues that differ from those in the wild-type substrate TLP are in bold type.  $^\dagger$ Lower limit of detection.

Many of the kinetic assays done by Rahul Kohli



Cyclization reactions: look for cyclization

Stock solutions: (3-5 mM)

3 mM ① FLPI (MW = 1484.6, TFA salt) 2.2 mg  $\rightarrow$  3 mM  
dissolve in 494  $\mu$ l H<sub>2</sub>O  $\Rightarrow$  3 mM

2.4 mM ② TLP3

2.43 mM stock solution previously prepared

3 mM ③ MIF (MW = 1354.4, 2 TFA salt) 3.5 mg  $\rightarrow$  3 mM  
dissolve in 862  $\mu$ l H<sub>2</sub>O  $\Rightarrow$  3 mM

3 mM ④ MIB (MW = 1109.1, 2 TFA salt) 3.2 mg  $\rightarrow$  3 mM  
dissolve in 962  $\mu$ l H<sub>2</sub>O  $\Rightarrow$  3 mM

1 mM ⑤ Phelac (MW = 504.7, 1 TFA salt) 1.4 mg  $\rightarrow$  1 mM  
dissolve in 930  $\mu$ l  $\Rightarrow$  1 mM

Enzymes:

	TE	prep	conc
①	TyCCTE	TE 89	40 $\mu$ M
②	PCP10TE	PCP1078	23 $\mu$ M
③	FenTE	FenTE 46	97 $\mu$ M

substrate: 2  $\mu$ M 10  $\mu$ M 50  $\mu$ M  
 enzyme: 50  $\mu$ M TE  
 200  $\mu$ M TE

97

### Cyclization conditions:

30 - 200  $\mu$ M peptide-SNAC  
 200  $\mu$ M TE  
 25 mM MOPS, pH 7.0  
 time = 2 min  
 } total volume of 400  $\mu$ L

Enzymes: PCPIOTE x 1  
 Fente x 1  
 TyccTE x 4

### Reactions

Rn	TE	Substrate	Amt. Peptide stock	Amt. H <sub>2</sub> O
1	TyccTE 200 $\mu$ M	TLP3 (320)	33 $\lambda$	287 $\lambda$
2	TyccTE 200 $\mu$ M	Phelac1 (320)	80 $\lambda$	240 $\lambda$
3	TyccTE 180 $\mu$ M	M17 400	27 $\lambda$	373 $\lambda$
4	TyccTE 180 $\mu$ M	M18 400	27 $\lambda$	373 $\lambda$
5	PCPIOTE 200 $\mu$ M	TLP3 (320)	33 $\lambda$	287 $\lambda$
6	Fente 180 $\mu$ M	FLP1 400	27 $\lambda$	373 $\lambda$

### Enzyme stocks:

10.6  $\lambda$  40  $\mu$ M TyccTE  
 189  $\lambda$  180  $\mu$ M TE buff.  
 200  $\lambda$  2  $\mu$ M TyccTE  
 5  $\lambda$  23  $\mu$ M PCPIOTE  
 45  $\lambda$  180  $\mu$ M TE buff.  
 50  $\lambda$  PCPIOTE, 2  $\mu$ M  
 4.3  $\lambda$  92  $\mu$ M Fente  
 136  $\lambda$  180  $\mu$ M TE buff.  
 200  $\lambda$  2  $\mu$ M Fente  
 TE buff = 10 mM MOPS, pH 7,  
 10 mM NaCl


### Reactions

40  $\lambda$  250 mM MOPS, pH 7  
 320  $\lambda$  H<sub>2</sub>O + peptide-SNAC  
 40  $\lambda$  TE stock  
 2 min at RT  
 - then + 25  $\lambda$  17% TFA/water quench  
 - freeze LN<sub>2</sub>  
 - + 80  $\lambda$  CH<sub>3</sub>CN (17% F/w)



↓ = new products

HPLC of Samples



# Reaction of $\beta$ -series w/ TycC TE.

## Stock soln of $\beta$ -series peptides.

	Amt	MW	DDW	[Final]
$\beta 6$	2.1 mg	1209.4	0.347 $\mu$ L	5 mM
$\beta 8$	1.6 mg	1403.6	0.456 $\mu$ L	2.5 mM
$\beta 12$	2.0 mg	1827.9	438 $\mu$ L	2.5 mM
$\beta 14$	1.6 mg	2026.0	790 $\mu$ L	1 mM

$\hookrightarrow$  purity estimated @ 50%, thus true final [A] near 500  $\mu$ M

Rxn:

Make 500  $\mu$ M stock of each peptide - SNAC.

stock =  
8.73 mM  $\leftarrow$

$\beta 6$	dilute	10 $\mu$ L $\rightarrow$	100 $\mu$ L
$\beta 8$	dilute	20 $\mu$ L $\rightarrow$	100 $\mu$ L
$\beta 10$	dilute	5.7 $\mu$ L $\rightarrow$	100 $\mu$ L
$\beta 12$	dilute	20 $\mu$ L $\rightarrow$	100 $\mu$ L
$\beta 14$	use	stock. (purity 50% $\rightarrow$	$\sim$ 500 $\mu$ M)

Rxn	(500 $\mu$ M) peptide	250 mM NaPS, pH 7	DDW	+ 2 $\mu$ M TycC TE to initiate	
$\beta 6$ TE	40 $\mu$ L	40 $\mu$ L	280 $\mu$ L	40 $\mu$ L $\rightarrow$	1' - Add 25 $\mu$ L 1.7% TFA Flash freeze
$\beta 6$ cont	"	"	"	-	
$\beta 8$ TE	"	"	"	40 $\mu$ L $\rightarrow$	1' - Add 25 $\mu$ L 1.7% TFA Flash freeze
$\beta 8$ cont	"	"	"	-	
$\beta 10$ TE	"	"	"	40 $\mu$ L $\rightarrow$	1' - Add 25 $\mu$ L 1.7% TFA FF
$\beta 10$ cont	"	"	"	-	
$\beta 12$ TE	"	"	"	40 $\mu$ L $\rightarrow$	1' - Add 25 $\mu$ L 1.7% TFA FF
$\beta 12$ cont	"	"	"	-	
$\beta 14$ TE	"	"	"	40 $\mu$ L $\rightarrow$	1' - Add 25 $\mu$ L 1.7% TFA FF
$\beta 14$ cont	"	"	"	-	

For all cont rxn, add first 25  $\mu$ L 1.7% TFA.

then add 40  $\mu$ L 2  $\mu$ M TycC TE  
Flash freeze instantly.

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**